KliniMark CLIP

REF: 105880900N

INSTRUCTIONS FOR USE

Manufacturer's contact details

KLINIKA Medical GmbH
Am Dornbusch 29, 61250 Usingen
Germany
Tel.: +49 6081 9139-0
info@klinika-medical.de

Feedback

In the event of serious incidents occurring in connection with this medical product, please report these to us as well as to the responsible authorities.

Feedback on the instructions for use

Your opinion is important to us. Feel free to share your suggestions and criticism about these instructions with us. We will assess your feedback and consider it for the next version of the instructions.

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
**Scope**

These instructions for use apply to the following products:

Article: KliniMark CLIP
Article number: REF: 105880900N

Important notice

Improper handling or use not in accordance with the intended purpose may result in hazards and damage. Please read and follow these instructions for use carefully. Keep these instructions for use readily accessible. Follow the safety instructions to prevent injury or device malfunction. If you have any questions about the contents of these instructions for use or the use of the product, please contact the KLINIKA team.

Product description

The KliniMark CLIP is a sterile, disposable product consisting of a puncture cannula (18G, 12 cm long) and a preloaded, non-absorbable Nitinol clip marker. The clip marker has an omega V shape with dimensions of 5.9 x 2.1 mm and a diameter of 4.5 mm. The cannula has 1 cm depth markings and a beveled tip for tissue insertion. The product is intended for long-term use.

Operating principle

The preloaded, non-magnetic Nitinol clip is introduced into the target area through an 18G cannula. The clip is deployed in the tissue by disengaging the spacer. Nitinol can be clearly seen on ultrasound and X-ray imaging. The cannula is inserted into the target area under imaging guidance, such as mammography or ultrasound. The marker is deployed by disengaging the spacer on the mandrel and advancing the mandrel.

The position of the marker is documented. The cannula is then withdrawn. Marking the tumor or other tissue areas enables subsequent localization or monitoring (e.g., in neoadjuvant therapy). The KliniMark clip's Omega-V shape enables reliable three-dimensional localization.

Information about the material

The Clip is made of Nitinol (a nickel-titanium alloy) and the cannula of steel.

Intended use

The KliniMark CLIP is a medical device intended for use in diagnosis, treatment, and monitoring. The KliniMark CLIP is used for tissue marking of breast tissue and lymph nodes. It may be used before or after chemotherapy or radiation therapy to mark the site of a biopsy, the location of a resected tumor, or areas with ambiguous or inconclusive findings.

Indications

The markers are indicated for use in lesions of the breast tissue and lymph nodes before or after chemotherapy or radiation therapy or at the site of a biopsy or a resected tumor. The KliniMark Clip can also be used for preoperative tissue marking in breast tissue and lymph nodes.

Contraindications

- Nitinol is a nickel-titanium alloy, so the clip may only be used in patients with a known severe nickel allergy if specifically approved by a physician.
- The marker must not be used in the heart, central nervous system or central circulatory system.
- The cannula is not MRI compatible and must not be deployed in an MRI environment.

MRI safety

Application system: The cannula is not MRI compatible, and must not be used in an MRI environment.

Clip marker: The KliniMark Clip has been tested for safety in magnetic resonance imaging (MRI) environments and has been found to comply with applicable ASTM standards. Based on the test results, the KliniMark Clip is considered MRI compatible under the following conditions and may be safely used in MRI devices:

Magnetic field strength: 1.5 Tesla and 3 Tesla
Maximum spatial field gradient: 20 T/m
Maximum whole-body specific absorption rate (SAR): ≤ 2 W/kg

Please note that the test results only apply to the above-mentioned conditions. Other MRI parameters may influence the behavior of the clip. Before any MRI examination, the patient's individual health status and any other implanted devices should be taken into account.

Clinical benefit

The clip's design ensures stable retention in tissue and significantly reduces the risk of unintended migration, supporting its intended use for localization.

Intended patient population

Adults. The product must not be used on children, pregnant women or breastfeeding women

Intended users

Medical staff

Intended use environment

Medical treatment room, operating room

Use

1. Before use, ensure that the packaging is intact and unopened. Please check the expiry date.
2. Open the packaging at the tear tab with the foil side facing up. Peel back the foil and remove the product using aseptic technique.
3. Disinfect the procedure site and, if necessary, administer local anesthesia.
4. Remove the needle cover.
5. Under mammographic or ultrasound guidance, puncture the target area with the cannula.
6. Deploy the marker by disengaging the spacer from the mandrel and advancing the mandrel.
7. The clip unscrews from the cannula holder on the side opposite the blue marker.
8. Document the position of the marker. Then remove the cannula.
9. Treat the wound.

10. After treatment, the cannula must be disposed of properly in a suitable cannula container.

Storage and shelf life

Do not use the product if the packaging is open, damaged, or wet. Store in a dry place, away from direct light and heat. Do not store or use the product at temperatures below +5°C. The product's sterility is assured if the packaging is intact. The KliniMark CLIP has a shelf life of three years when stored in clean, cool, dry, and dark conditions.

Warnings


- The KliniMark CLIP should only be used by qualified physicians with sufficient experience and expertise in percutaneous soft tissue marking.
- This product is intended for single use only. It must not be resterilized or reprocessed after a single use. Resterilization or reprocessing may degrade the device materials or compromise its safety. Resterilization may lead to serious infections, cross-contamination, and injury to the patient. The preloaded Nitinol marker in the cannula is non-absorbable and intended for single use only. KLINIKA Medical GmbH disclaims all liability for resterilization or reuse of any component of the product.
- The marker must not be used in the heart, central nervous system or central circulatory system.
- Only the CLIP itself may remain in the tissue.
- The cannula is not MRI compatible and must not be deployed in an MRI environment.
- Please only use the KliniMark CLIP before the expiry date and only if the packaging is intact and unopened. Failure to do so will compromise the sterility of the product.

Potential adverse effects

None known

Residual risks

- Damage to the clip surface may result in a nickel-related allergic reaction.
- In rare cases, the clip may migrate within the tissue and thus fail to mark the target area.
- The clip may deform, which can reduce its visibility on ultrasound.
- The cannula may become hot during MRI and cause burns.
- Improper deployment may result in injury, hematoma, or pain.
- Bleeding may occur if a blood vessel is damaged by the cannula or clip.
- Reuse may lead to infections and cross-contamination.

SSCP

A Summary of Safety and Clinical Performance (SSCP) for this implantable device, prepared in accordance with Article 32 of EU Regulation 2017/745, is available at www.klinika-medical.de.

Product information

Design	KliniMark CLIP, sterile steel cannula with Nitinol marker
Diameter (gauge)	18
Length of cannula (mm)	120
Cannula material	Steel
Size of Clip (mm)	4.9 x 2.1
Diameter of Clip (mm)	4.5
Clip material	Nitinol
PU	2 x 5
Article no.	105880900N
Sterilization method	EO / EtO

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