FR

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Instructions for use of bottles, nipples, breast shields and accessories sterilised with ethylene oxide (References: NCB; NCT)

1. Indication:

The bottles, nipples, breast shields and accessories are single-use devices, for medical use only. They are indicated in the collection, storage and administration of breast milk to premature infants or infants suffering from serious illnesses.

Set of breast shields:

- NCB4....: Milk collection with Kitett breast pump
- NCB5....: Milk collection with Spectra breast pump
- NCB6....: Milk collection with Medela or Ameda breast pump

NOESY...
Set of *nipples*:

- NCT2.../NCT3.../NCT8...: To be screwed onto a bottle for administration of breast milk.

Accessories:

- NCT010A: Milk cup.

Users: These devices are intended to be used by healthcare professionals in a hospital setting and at home for mothers who donate their milk.

Target population: NCB feeding bottles... - Neonatology (>2kg); NCT nipples and accessories... - Neonatology (>2kg); Set of breast shields NCB4.../NCB5.../NCB6... - Adult (>58kg).

2. Composition:

- NCB4.../NCB5.../NCB6...: PP breast shield; PVC tube; PVC cup; PP bottle (if there is one).
- NCT2.../NCT3.../NCT8...: SEBS nipples; SEBS ring; PP cover.
- NCT010A: PP milk cup.

PE= Polyethylene / HDPE= High-density polyethylene / PP= Polypropylene / SEBS= Styrene-ethylene-butylene-styrene / PVC = polyvinyl chloride

3. Performance:

- Devices sterilised with ethylene oxide to prevent risks of infection.
- Sealed devices for the collection, storage and administration of breast milk.

4. Instructions for use:

Prerequisites before use

- Check the integrity and cleanliness of the packaging and the device before use; if there are any defects, do not use the device.
- Verify all connections and assemblies before use.
- Check the expiry date before use.
- Carefully open the cavity, taking care with the peel-off direction.

Instructions for use

Please refer to the recommendations and hygiene procedures in force in the facility for all steps of use.

Collection (Ref.: NCB4.../NCB5.../NCB6...)

- Screw the bottle on the breast shield, connect the tube to the breast pump or vacuum. Close the cap of the breast shield to enable milk collection. Place the breast in the breast shield to enable milk collection.

Administration (Ref.: NCT2.../NCT3.../NCT8.../NCT010A)

- Rewarm the milk in the dry bottle warmer or using the water bath, according to hospital protocol.
- Screw on the nipple ring (previously assembled with the nipple) or feeding cap onto the bottle to initiate administration of the milk.
- Stir the milk once heated to ensure even heat distribution and check the temperature before initiating feeding.

Duration of use

The devices should be used during the time of care, according to hospital protocol.

5. Precautions for use and contraindications

- After use, the devices cannot be reused and must be disposed of in accordance with the procedures in force.

Contraindications

- The 60-ml bottles cannot be pasteurised or frozen.
- Do not mix freshly collected milk with frozen milk.
- Do not break the cold chain. Do not refreeze thawed milk.
- Do not let the milk rewarm at room temperature due to the risk of microbial development.
- Use of a microwave to rewarm the contents of the bottle is prohibited.
- Always stir the food once heated to ensure even heat distribution and check the temperature before serving.

6. Risks and adverse events:

- There are no known adverse effects related to the use of these medical devices.
- In the event of a serious incident with the medical device, the user must inform the manufacturer and the relevant health authority of the country where the incident occurred.

The medical devices involved must absolutely be kept so that they can be examined by an expert.

7. Storage - Handling - Disposal:

Storage: The device must be stored in its original packaging, in a clean, dry area that is away from the light and at a temperature of between +5 and +40°C.


Shelf life: Refer to the expiry date that appears on the label.

Waste disposal: To dispose of this device, please refer to the regulation applicable to the disposal of hazardous waste, biological waste or infectious medical waste

8. Single use

CAIR LGL did not design this device to be reprocessed, reused or resterilised and, therefore, cannot certify that reprocessing (e.g.: cleaning or sterilisation) will maintain the device's structural integrity and guarantee the safety of the patient and/or user.

9. Glossary of symbols

 Do not reuse. Intended only for single use.	 Do not use if the packaging is damaged.	 Expiry date.	 Always read the instructions for use.	 Read the instructions for use relating to warning information.	 Reference number.	 Batch code.	 Manufacturer	 Manufacturing date.
 Medical Device	 Do not resterilise	 Peel-off direction	 Unique sterile barrier system	 Sterilised with Ethylene Oxide	 Store away from sunlight	 Fears humidity	 Quantity	 Temperature limit (5°C – 40°C)
 Not made with latex		 Not made with DEHP		 Not made with Bisphenol A		 Not made with PVC		 Internal diameter

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