KliniMark CLIP PRO

REF: 205880900C

INSTRUCTIONS FOR USE

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Feedback

In the event of serious incidents occurring in connection with this medical product, please report these to us as well as to the responsible authorities.

Feedback on the instructions for use

Your opinion is important to us. Feel free to share your suggestions and criticism about these instructions with us. We will assess your feedback and consider it for the next version of the instructions.

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Scope

These instructions for use apply to the following products:

Article: KliniMark CLIP PRO
Article number: 205880900C

Important notice

Improper handling and use can result in hazards and damage. We therefore ask you to read through these instructions for use and follow them exactly. Keep them close at hand at all times. To avoid personal injury and damage to property, please also observe the safety instructions. If you have any questions about the contents of these instructions for use or the use of the product, please contact us. Your KLINIKA team

Product description

The KliniMark CLIP PRO is a sterile, disposable product consisting of a puncture cannula (18G, 15,5 cm long) and a preloaded, non-absorbable Nitinol clip marker. The clip marker has an omega shape with dimensions of 5.9 x 2.1 mm and a diameter of 4.5 mm. The cannula has space markings (1 cm) for measuring the penetration depth and a bevelled tip for insertion into the tissue. The product is intended for long-term use.

Features

Using imaging feedback, such as mammography or ultrasound, the cannula is guided to puncture the target area. This ensures the exact placement of the marker. The applicator slider is first moved to the loading position. This prepares the clip inside the cannula without prematurely discharging it. When the slider is advanced more, the Nitinol clip is gently applied to the tissue. This can be easily done with a single hand, ensuring precise control during the procedure.

The Nitinol clip is made of non-magnetic material and is clearly visible on both ultrasound and X-ray images. This allows the position of the marker to be checked at any time after placement. The position of the marker is documented. Then the cannula is removed. It is helpful to mark tumours or other areas in this way so that they can be relocated or observed

(e.g. in neoadjuvant therapy). Due to its shape (omega V), the KliniMark CLIP PRO is easy to relocate three-dimensionally.

Information about the material

The clip is made from Nitinol (nickel-titanium alloy), the cannula from steel.

Intended use

The KliniMark CLIP PRO is a medical device used for diagnosis, treatment and monitoring in medical procedures. KliniMark CLIP PRO is used for tissue marking in breast tissue and lymph nodes: Before or after chemotherapy or radiotherapy to mark the location of a biopsy removal or a removed tumour, in the case of unclear findings.

Indication

The markers are used in lesions of the breast tissue and lymph nodes before or after chemotherapy or radiotherapy or at the site of a biopsy or a removed tumour. KliniMark Clip PRO can also be used for preoperative tissue marking in breast tissue and lymph nodes.

Contra-indications

- Nitinol is a nickel-titanium alloy and should therefore, in the case of a known, pronounced nickel allergy, only be used after a doctor's recommendation.
- The marker must not be used on the heart, central nervous system or central circulatory system.
- The cannula is not MRI capable and must not be placed under MRI.

MRI safety

Application system: The cannula is not MRI safe and must not be used in an MRI environment.

Clip marking: The KliniMark CLIP PRO has been tested for safety in magnetic resonance imaging (MRI) environments and evaluated in accordance with ASTM standards. Based on the test results, the KliniMark CLIP PRO is considered MRI-conditional safe under the following conditions and is safe to use in MRI machines:

Magnetic field strength: 1.5 Tesla and 3 Tesla
Maximum spatial field gradient: 20 T/m
Maximum whole-body specific absorption rate (SAR): ≤ 2 W/kg

Please note that the test results are only valid for the conditions mentioned above. Other MRI parameters may affect the behavior of the clip. Before any MRI examination, the patient's individual health conditions and any other implants that may be present should be taken into account.

Clinical benefit

The clip's specific shape ensures a secure grip within the tissue and significantly reduces the risk of unwanted clip migration, thus increasing the clinical benefit for the patient.

Designated patient group

Adults. The product must not be used on children, pregnant women or breastfeeding women.

Designated users

Medical staff

Designated environment for use

Medical treatment room, surgery

Use:

1. Before use, it must be checked that the packaging is closed and not damaged. Please check the expiry date.
2. Open the packaging at the tear-off flap with the foil side facing up. Remove the foil and take out the product aseptically.
3. Disinfect the area of the procedure and, if necessary, perform local anaesthesia.
4. Remove the protective cap from the needle.
5. Under an imaging procedure, such as mammography or ultrasound, puncture the target area with the cannula.

6. Push the applicator slider to the release position. This prepares the clip without discharging it. Check that the slider is properly engaged before proceeding with the application.
7. Push the slider forward in a controlled manner to apply the Nitinol clip to the tissue. One-handed operation ensures precise control.
8. Document the position of the marker. Then remove the cannula.
9. Treat the wound.
10. After treatment, the cannula must be disposed of properly in a suitable cannula container.

Storage and storage life

The product must not be used if the packaging is open, damaged or damp. Store in a dry place, away from direct light and heat.

Do not store or use the product at temperatures below +5 °C. Sterility is guaranteed by the complete packaging. The KliniMark CLIP PRO has a storage life of three years in clean, cool, dry and dark conditions.

Warning notes

- The KliniMark CLIP PRO should only be used by qualified physicians who have sufficient experience and expertise in the percutaneous marking of soft tissue.
- This product was manufactured expressly for single use. It must not be resterilised or recycled after a single application. Resterilisation or recycling may result in material defects of the product or the medical device becoming no longer safe for use. Resterilisation may lead to serious infections, cross-contamination and injury to the patient. The preloaded Nitinol marker in the cannula is non-absorbable and intended for single use only.

KLINIKA Medical GmbH rejects any responsibility in the case of resterilisation and reuse of even one of the product components.

- The marker must not be used on the heart, central nervous system or central circulatory system.
- Only the clip may remain in the tissue.
- The cannula is not MRI capable and must not be placed under MRI.
- Please only use the KliniMark CLIP PRO before the expiry date and with undamaged, unopened packaging. Otherwise the sterility of the product cannot be guaranteed.

- If the product is not applied properly, this can result in injury, hematomas and pain for the patient.
- Bleeding may occur if a blood vessel is damaged by the cannula or clip.
- Reuse may lead to infections and cross-contamination.

SSCP

A summary report of the safety and clinical performance pursuant to EU 2017/745, Art. 32 on implantable products (SSCP) will be made available on the website www.klinika-medical.de.

Side effects

None known

Residual risks

- If the surface of the clip becomes damaged, this may lead to allergic reactions (nickel).
- In rare cases, the clip may migrate in the tissue and thus fail to mark the desired area.
- Crosslinking of the clip may occur which can reduce visibility in the ultrasound image.
- The cannula may become hot in MRI and cause burns.

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